

Suspend

Glider recliner



Featuring a recliner that is a glider in the upright position. The gliding motion automatically stops as the chair starts to recline. The footrest and first recline position are activated by the lever, push on the arms to reach the fully reclined position.



Features and Performance



Design Features

1. Continuous upholstery on foot rest (no gaps) for comfort, safety and ease of cleaning
2. Composite foam seat with dynamic suspension system for long term comfort
3. Polyester or wood arm caps resistant to hospital cleaners
4. Glider in the upright position. The gliding motion automatically stops as the chair starts to recline.

Performance

- 350 lb capacity
- Heavy duty mechanism completely covered by all inclusive 12 year warranty

Environmental Sustainability

- HHI Compliant
- Air Quality Certification BIFMA 207.7HI

Back Styles



straight back 628-S



flared back 628-F

*Images shown may not reflect style

Options and Accessories



Profile tablet

The profile tablet provides the patient with a generous and sturdy multi-purpose surface. When not in use it easily folds flush against the side of the chair.



Folding side tablet

A molded plastic tablet with raised edges that folds flush on the side of the arm panel.



Swiveling side tablet

Swiveling table that can be positioned in any position then locked. The swiveling table is also removable and available in 3 sizes.



Folding arm

Arm panel that can be folded to access the patient from the side of the chair. This is particularly useful for transferring patients to and from the bed.



Softgrip push bar

A tubular push bar with a softgrip cover for pushing the recliner.



Ergotron arm

Supports TV/DVD players and may be adjusted into a multitude of positions.



Drainage bag holder

A retractable loop for foley bags. It can be placed anywhere on the arm panel of the recliner.






Heated seat and lumbar

The heated seat and lumbar provides therapeutic warmth and eliminates the need of warming blankets. The optional massaging action provides three variations: wave, pulsating, and vibrating.

Dimensions

Glider recliner

Model	Short Description	W	D	RD	H	SW	SD	SH	AH	WT	YDG	WT CAP
628s-35	recliner, glider recliner straight back	32	40	66	41	24	22	18	25	160 .lbs	7	350 .lbs
628f-35	recliner, glider recliner flared back	32	40	66	41	24	22	18	25	160 .lbs	7	350 .lbs



*Images shown may not reflect style