Suspend

Motorized bariatric recliner with patient lift assist



The patient lift recliner is designed to enable better, safer interactions between caregivers and patients. With a 750 lb. capacity and motorized movements, the bariatric lift chair fosters a feeling of independence and will assist larger patients to their feet.





Features and Performance





Design Features

- Continuous upholstery on foot rest (no gaps) for comfort, safety and ease of cleaning
- 2. Composite foam seat with dynamic suspension system for long term comfort
- 3. Shielded mechanism for safety
- 4. Polyester or wood arm caps resistant to hospital cleaners
- 5. At the touch of a button, the patient may recline to a comfortable position
- 6. Assist larger patients to their feet with a touch of a button

Performance

- 750 lb capacity
- Heavy duty mechanism completely covered by all inclusive 12 year warranty

Environmental Sustainability

- HHI Compliant
- Air Quality Certification BIFMA 207.7HI

Back Styles



straight back 628-S



flared back 628-F



Options and Accessories



Profile tablet

The profile tablet provides the patient with a generous and sturdy multi-purpose surface. When not in use it easily folds flush against the side of the chair.



Folding side tablet

A molded plastic tablet with raised edges that folds flush on the side of the arm panel.



Swiveling side tablet

Swiveling table that can positioned in any position then locked. The swiveling table is also removable and available in 3 sizes.



Folding arm

Arm panel that can be folded to access the patient from the side of the chair. This is particularly useful for transferring patients to and from the bed.



Softgrip push bar

A tubular push bar with a softgrip cover for pushing the recliner.



Ergotron arm

Supports TV/DVD players and may be adjusted into a multitude of positions.



Drainage bag holder


A retractable loop for foley bags. It can be placed anywhere on the arm panel of the recliner.





Heated seat and lumbar

The heated seat and lumbar provides therapeutic warmth and eliminates the need of warming blankets. The optional massaging action provides three variations: wave, pulsating, and vibrating.

Dimensions







Motorized bariatric recliner with patient lift assist

Model	Short Description	W	D	RD	Н	SW	SD	SH	АН	WT	YDG	WT CAP
628s-09m-750	recliner, Bariatric lift chair	31	37	76	46	26	20	20	26	264 lbs	9	750 lbs
628f-09m-750	recliner, Bariatric lift chair	31	37	76	46	26	20	20	26	264 lbs	9	750 lbs

